



Spot Vital Signs

Product Specifications

Dimensions: h: 9.70, l: 5.72", w: 4.73"

h: 24.64 cm, l: 14.53 cm, w: 12.01 cm

Weight: 4.5 lbs. (2.04 kg)

Blood Pressure Accuracy: Meets or exceeds SP10-1992 AAMI standards

Cuff Pressure Range: 0-300 mmHg

Systolic Determination: Max: 250 mmHg Min: 60 mmHg

Diastolic Determination: Max: 160 mmHg Min: 30mmHg

Pulse Rate Accuracy:

SpO2 Module Heart Rate: ± 3.0 bpm

Blood Pressure Algorithm: ± 5.0 bpm

Pulse Rate Determination: Max: 245 bpm Min: 25 bpm

Temperature Accuracy: Meets or exceeds ASTM E112-86 standards

Pulse Oximeter Accuracy:

Nellcor Puritan Bennett Pulse Oximeter Module: $\pm 3\%$ in the range of 70-100% Oxygen saturation (1 standard deviation) < 70% unspecified by the OEM

Power Requirements: North American Version: 120 VAC, 60 Hz

International Version: 230-240 VAC, 50 Hz

Internal Battery: Lead Acid, with external charger

A complete listing of specifications is available on request.