



Meena Medical Equipment, Inc. • 1904 Industrial Blvd., Suite 105 • Colleville, TX 76034 • 1-888-225-2502

Kendall SCD System 6325 Manufacturer Specifications

MODEL:	SCD System 6325
PUMP SPECIFICATIONS:	
Press range, mm hg (PUMP SPECIFICATIONS):	45
CYCLE DURATION, SEC (APERTURE):	71
Inflation (PUMP SPECIFICATIONS):	11
Deflation (PUMP SPECIFICATIONS):	60
LINE POWER, VAC (TEMPERATURE):	115
Power, w:	100
Power switch (TUBING):	Yes
USES (STAINING POSTS):	DVT prophylaxis
LINING:	TED recommended **
DISPLAY:	LED
Alert (DISPLAYS):	Audible and visual - auto decompresses with pressure fault
TUBING:	
Connector type (PROBE):	Snap fit
Correct alignment:	
Indicator (ACCESS):	Alignment of arrows/ keyed
PNEUMATIC SLEEVE:	
Form (PNEUMATIC SLEEVE):	Thigh - full-leg SM,REG, XL; knee length
MATERIAL:	PVC or polyolefin with absorbent lining
Closure (PNEUMATIC SLEEVE):	Velcro
TYPE/NO. OF GARMENTS (TUBING):	3-chamber gradient sequential/2
Sleeve Cooling Function	Yes
PUMP SIZE, HXWXD CM (IN) (PUMP SPECIFICATIONS):	15.2x 26.3 x 11.4 (6x 10.34 x 4.5)
WEIGHT, kg (lb) (DISPLAY):	4(8.9)
OTHER ATTRIBUTES:	Unit has a cooling
Bed Hook	Yes
Power Cord Storage	Yes
FDA CLEARANCE:	Yes



**Meena Medical
Equipment, Inc.**