



Meena Medical Equipment, Inc. • 1904 Industrial Blvd., Suite 105 • Colleville, TX 76034 • 1-888-225-2502

Kendall SCD System 7325 Manufacturer Specifications

MODEL:	7325 SCD Response
PUMP SPECIFICATIONS:	
Deflation (PUMP SPECIFICATIONS):	60, variable; vascular refill detection
Inflation (PUMP SPECIFICATIONS):	11
Pressure range, mm Hg (PUMP SPECIFICATIONS):	45 ankle, 40 calf, 30 thigh
CYCLE DURATION, SEC (APERTURE):	71variable
PUMP SIZE, HXWXD CM (IN) (PUMP SPECIFICATIONS):	15.2 x 26.3 x 11.4 (6 x 10.3 x 4.5)
WEIGHT, kg (lb) (DISPLAY):	4 (8.9)
USES (STAINING POSTS):	DVT prophylaxis
DISPLAY:	LED
Alert (DISPLAYS):	Audible and visual - auto decompresses with pressure fault
LINE POWER, VAC (TEMPERATURE):	115/230
Power switch (TUBING):	Yes
Power, w:	100
TUBING:	
Alignment signal (TUBING):	Alignment of arrows/keyed
Connector type (PROBE):	Snap on
Non-IV-compatible connector (TUBING):	Yes
PNEUMATIC SLEEVE:	
Closure (PNEUMATIC SLEEVE):	Hook and loop
Form (PNEUMATIC SLEEVE):	thigh (S, M, L), knee (M, L, XL)
MATERIAL:	PVC with nonwoven liner
TYPE/NO. OF GARMENTS (TUBING):	3-chamber gradient sequential/2
Bed Hook	Yes
Power Cord Storage	Yes
OTHER ATTRIBUTES (Interference compensation):	Vascular Refill Detection technology determines venous refill time of individual patients and sets deflation time accordingly; passive cooling system of air pumped through sleeve. Meets requirements of UL 2601.
FDA CLEARANCE:	Yes



**Meena Medical
Equipment, Inc.**