



**Meena Medical
Equipment. Inc.**

Kendall SCD Express 9525 Manufacturer Specifications

MODEL:	9525 SCD Express
PUMP SPECIFICATIONS:	
CYCLE DURATION, SEC (APERTURE):	71variable
Pressure range, mm Hg:	45 ankle, 40 calf, 30 thigh
Inflation (PUMP SPECIFICATIONS):	11
Deflation (PUMP SPECIFICATIONS):	60, variable; vascular refill detection
PUMP SIZE, HXWXD CM (IN):	15.7 x 17.7 x 11.4 (6.2 x 7 x 4.5)
WEIGHT, kg (lb) (DISPLAY):	1.6 (3.5)
USES (STAINING POSTS):	DVT prophylaxis
Alert (DISPLAYS):	Audible & visual - auto decompresses with pressure fault
DISPLAY:	LED
LINE POWER, VAC (TEMPERATURE):	115/230; also lithium ion battery
Power switch (TUBING):	Yes
Power, w:	100
PNEUMATIC SLEEVE:	
Closure (PNEUMATIC SLEEVE):	Hook and loop
Form (PNEUMATIC SLEEVE):	Thigh (S, M, L), knee (M, L, XL), foot (M, L)
MATERIAL:	PVC with nonwoven liner
TUBING:	
Alignment signal (TUBING):	Alignment of arrows/keyed
Connector type (PROBE):	Snap on
Non-IV-compatible connector (TUBING):	Yes
TYPE/NO. OF GARMENTS (TUBING):	3-chamber gradient sequential/2
OTHER ATTRIBUTES:	Vascular Refill Detection technology determines venous refill time of individual patients and sets deflation and cycle time accordingly; all-in-one controller allows both leg and foot compression; SCD Express Kambia sleeve allows users to tear away
FDA CLEARANCE:	Yes
MARKETING REGION:	Worldwide