

Model FP3B Handheld Doppler



- **The Gold Standard**

A classic, rugged design recognized for reliability and consistency

- **Practical**

Easy to use configuration incorporating the Doppler and probe into one device

- **Confidential**

Using a standard stethoscope allows discreet screening of fetal condition without alarming the patient

- **Unique**




Narrow beam Doppler signal allows for increased sensitivity

Model FP3B **Fetal Doppler**



Model FP3B

Fetal Doppler Ultrasound Stethoscope®

ORDER INFORMATION

	Description	Order Number	Frequency	Contents - Applications /Specifications
	Model FP3B Fetal Doppler	101-0153-010	2.25 MHz	Includes: Doppler, stethoscope headset, gel, battery, soft carrying case. Applications: Detect fetal life early in pregnancy. Assess the rate and rhythm of the fetal heart.
	FetalPulse Plus™ Fetal Doppler System	101-0257-010	2.25 MHz	Includes: Doppler, stethoscope headset, FetalCalc™ speaker/heart rate display, gel, batteries, hard carrying case. Applications: Detect fetal life early in pregnancy. Assess the rate and rhythm of the fetal heart.
	Model FP3B Fetal Doppler & SA2 Pocket Speaker	101-0257-020	2.25 MHz	Includes: Doppler, stethoscope headset, SA2 Pocket Speaker, gel, batteries, soft carrying case. Applications: Detect fetal life early in pregnancy. Assess the rate and rhythm of the fetal heart.

OPTIONAL DOPPLER ACCESSORIES

	Model SA2 Pocket Speaker	101-0052-020	n/a	Includes: Speaker, cord to connect to Doppler. Specifications Power Source: One 9V alkaline battery Size: 4.2 x 2.5 x 1.1 in (10.8 x 6.3 x 2.7 cm) Weight: 6 oz (170 gm)
	Model SA6 Tabletop Speaker	101-0240-010	n/a	Includes: Speaker, cord to connect to Doppler. Specifications Power Source: Four C-cell alkaline batteries Size: 5.5 x 3.5 x 6.5 in (15 x 8.9 x 16.5 cm) Weight: 2.25 lb (1 kg)

FetalCalc™ SPECIFICATIONS

FetalCalc is a component of FetalPulse Plus. It is also available as an add-on accessory (Order Number 101-0238-010) to the Model FP3B.

Outputs	One 3.5 mm jack to connect to Doppler
Controls	Automatic ON/OFF, speaker disable switch
Power Source	One standard 9V alkaline battery
Weight	6.5 oz. (178 gm)
Dimensions	5.75 x 2.75 x 2 in (14.5 x 7 x 2.5 cm)
Construction	Rugged ABS plastic case, solid state circuitry

Additional Accessories: Replacement cords, stethoscopes and repair kits are also available. Contact your Sales Representative for part numbers and pricing.

Model FP3B SPECIFICATIONS

Ultrasonic Frequency	2.25 MHz
Outputs	Two 2.5 mm jacks to drive one or two headsets, one headset and a speaker or recording device.
Controls	Volume control knob and push button ON/OFF
Power Source	One standard 9V alkaline battery
Weight	4.7 oz. (133 gm)
Dimensions	4.95 x 2.5 x 1.0 in (12.6 x 6.4 x 2.5 cm)
Construction	Rugged ABS plastic case, integral transducers, solid state circuitry
Service:	CooperSurgical, Inc. maintains a full service facility which can promptly repair products returned to the factory.
Warranty:	MedaSonics Dopplers are covered by a 5-year warranty on parts; 1-year warranty on labor. <i>Call Customer Service for details.</i>

MEDASONICS®

a division of
COOPERSURGICAL

95 Corporate Drive, Trumbull, CT 06611
800-243-2974 • 203-601-5200
Fax: 800-262-0105
www.coopersurgical.com