



**Meena Medical
Equipment, Inc.**

Stackhouse, AIRSAFE ST-2500 Manufacturer Specifications

MODEL	AIRSAFE ST-2500
FILTRATION	
Prefilter needed (FILTRATION)	Yes
Filter(s) (FILTRATION)	ULPA/charcoal
Filter life (FILTRATION)	100 hr (ULPA), 6 or 12 hr (charcoal)
Filter life signal (FILTRATION)	Visual
POWER NEEDED	120/240 V
DIMENSIONS (HXWXD) CM, (IN) (DISPLAY)	87.6x 41.9 x 50.8
WEIGHT, kg (lb) (DISPLAY)	31.8(70)
Type (CONFIGURATION)	Stand-alone mobile system
EVACUATION WAND (TUBING)	Yes
Esu/laser (TUBING)	ESU
Add-on/integrated (TUBING)	Yes
Tube diam, mm (in) (TUBING)	6(1/4)
Extendable (TUBING)	Yes
DISPLAY	LED
DATA DISPLAYED (REGULATORS)	Filter status,
FLOW LIMITS	
Setting, cfm (FLOW LIMITS)	50 @ 22 mm
Variable flow control (FLOW LIMITS)	Yes
SOUND LEVEL, DBA (DIMENSIONS (LXWXH) CM, (IN))	60
Sizes, mm (in) (TUBING)	6(1/4),22(7/8),
SAFETY MEASURES	Filters w/port caps
Footswitch (US GENERATOR)	Pneumatic
SPECIFIED USE (TUBING)	Removal of surgical smoke / fluid in OR setting
ADDITIONAL ATTRIBUTES (MICROMANIPULATORS)	Fluid canister